

HORIBA

Explore the future

Sulfur-In-Oil-Analyzer



SLFA-2800/2100



SLFA-20




Specifications

	SLFA-2100	SLFA-2800	SLFA-20
Measurement Principle	Fluorescence X-ray Analysis Method		
Samples Measured	Sulfur in petroleum products such as heavy oil, naphtha and crude oil.		
Measurement Range	0 to 9.999%		0 to 5%
Repeat Accuracy	5 ppm or less (with a 1% sample), 1.6 ppm or less (with a 0% sample)		15 ppm or less (with a 1% sample)
C/H Error	±50 ppm or less per C/H (with a 1% sample)		
Lower Detection Limit	5 ppm or less		20 ppm or less (3 times the standard deviation)
No. of Analytical Curves	5 (both automatic selection and manual setting are possible)		
Analytical Curve Order	Primary or Secondary (both automatic selection and manual setting are possible)		
Calibration	Calibration made using any standard sample (2 to 20 calibration points)		
Sample Cell	Dedicated Liquid Sample Cell		
Sample Weight	4 to 10 ml		
Measurement Time	User-set from 10 to 600 seconds		10, 30, 100, 300, 600 seconds
Max. No. of Samples	1	8	1
Ambient Temperature	5 to 35°C		5 to 40°C
Ambient Humidity	30% to 85% (Relative Humidity)		80% maximum relative humidity at temperatures from 5°C to 31°C. Linear decrease to 50% relative humidity at temperatures from 31°C to 40°C.
Spectrum Measurement	Includes spectrum measurement functions, Energy Axis: 0-10 keV, Strength Value Axis: set automatically		
External Output	RS-232C output, output during measurement, automatic calibration, and spectrum measurement		
Electric Power Supply	AC 100 to 240 V ±10% 50/60 Hz		
Electric Power Consumption	150 VA		80 VA
External Dimensions	430(W) × 500(D) × 230(H) mm 16.9(W) × 19.7(D) × 9.1(H) in		250(W) × 407(D) × 137(H) mm 9.8(W) × 16(D) × 5.4(H) in
Weight	approx 21 kg	approx 24 kg	8 kg

* Standards Applied to Device Measurement Functions: JIS K2541, JIS B7995 (Japan), ASTM D429 (USA), ISO 8754

Dimensional Outlines Unit: mm(in)



Horiba continues contributing to the preservation of the global environment through analysis and measuring technology.



⚠ Please read the operation manual before using this product to assure safe and proper handling of the product.

- The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.



HORIBA <http://www.horiba.com>

- **HORIBA, Ltd.**
Head Office: Miyahogashi, Kisshoin Minami-ku, Kyoto, Japan
Phone: 81 (75) 313-8123
Fax: 81 (75) 321-5725
- **Tokyo Sales Office**
1-7-8 Higashi-Kanda Chiyoda-ku, Tokyo, Japan
Phone: 81 (3) 3861-8231
Fax: 81 (3) 3861-8259
- **Beijing Representative Office**
Room No. 410, No. 33 Cheng Fang Street, Xicheng District, Beijing P.R.O.C. 100032
T.O.D. Beijing-9503-316
Phone: 86 10-66077630/5
Fax: 86 10-66077554
- **Taiwan Representative Office**
4F, No.23-3, Lane 525, Sec. 1, Kuang Fu Rd., Hsin-Chu, Taiwan, R.O.C.
Phone: (886) 3-5799143
Fax: (886) 3-5799164
- **HORIBA KOREA Ltd.**
112-6 Sogong-Dong Choong-ku, Seoul, Korea
Phone: 82 (2) 753-7911
Fax: 82 (2) 756-4972
- **HORIBA INSTRUMENTS Pte. LTD.**
31 Kaki Bukit Road 3, #06-12 Techlink Singapore 417818
Phone: 65 745-8300
Fax: 65 745-8155
- **HORIBA INSTRUMENTS INCORPORATED**
Irvine Facility: 17671 Armstrong Avenue Irvine, CA 92614, U.S.A.
Phone: 1 (949) 250-4811
Fax: 1 (949) 250-0924
- **Ann Arbor Facility**
5900 Hines Drive Ann Arbor, MI 48108 U.S.A.
Phone: 1 (734) 213-6555
Fax: 1 (734) 213-6525
- **HORIBA / STEC INCORPORATED**
1080 E. Duane, Suite. A Sunnyvale, CA 94086 U.S.A.
Phone: 1 (408) 730-4772
Fax: 1 (408) 730-8975
- **HORIBA GmbH**
Kaplanstrasse 5 A-3430 Tulln, Austria
Phone: 43 (2272) 65225
Fax: 43 (2272) 65230
- **HORIBA CZECHIA**
Organizacni slozka Praha Petrohradská 13 CZ-101 00 Praha 10, Czech Republic
Phone: 420 (2) 717-464-80
Fax: 420 (2) 717-470-64
- **HORIBA EUROPE GmbH**
Head Office: Hauptstrasse 108 65843 Sulzbach Germany
Phone: 49 (6196) 6718-0
Fax: 49 (6196) 641198
- **Leichlingen Facility**
Julius-kronenberg Strasse D-42799 Leichlingen Germany
Phone: 49 (2175) 8978-0
Fax: 49 (2175) 8978-50
- **HORIBA FRANCE**
Rue L. et A. Lumière Technoparc F-01630 St-Genis-Pouilly France
Phone: 33 (4) 50-42-27-63
Fax: 33 (4) 50-42-07-74
- **HORIBA SWEDEN**
Hertig CarlsVag 55-57 S-15138 Södertälje Sweden
Phone: 46 (8) 550-80701
Fax: 46 (8) 550-80567
- **HORIBA INSTRUMENTS LIMITED**
Kyoto Close Summerhouse Road Moulton Park, Northampton NN3 6FL, U.K.
Phone: 44 (1604) 542500
Fax: 44 (1604) 542699

Bulletin:HRE-00000

Printed in Japan XX-X(XX)00

Our new sulfur-in-oil analyzer combines the established functions with new functions to improve the ease of use and enlarge the application.

Today's world is facing many environmental problems. In order to prevent harmful pollution high precision analysis technology with a high degree of sensitivity, precision and repeatability gains importance. The effort is to decrease the sulfur content in diesel fuel, light and heavy oil. This is due to present regulations concerning the sulfur content, which become even stricter in the future coupled with the rising environmental consciousness.

HORIBA's SLFA-2800/2100 meets all demands for on innovative analytical technology.

New features that enlarge the application and make the use of the analyzer easier.

The SLFA-2800/2100 includes several convenient new features without increasing its size.

The software has been updated, measure range has been enhanced and the new sulfur-in-oil analyzer presents itself in the new Horiba design.

■ Software

The new software does not only simplify the use of the analyzer; it also extends its application with new functions.

Measurement units

For the measurements it is possible to choose between the units % and ppm.

Sample identification

The new SLFA-2800/2100 makes the identification of the samples very easy. Sample IDs can include numbers (0-9) and alphabetic characters (A-Z) up to ten digits.

Calibration

In this new sulfur-in-oil analyzer the calibration graphs have been extended. The analyzer selects automatically out of seven scales (0.002 / 0.01 / 0.05 / ... / 10 %). The graphs are available in both units, % and ppm.

The curve format, linear or parabolic, can be either picked manually or set automatically. It is possible to store up to five calibration curves. You can choose if the analyzer should select the most suitable curve for you or if you want to do it manually.

■ Measure range

This new model can measure a wider range of sulfur content. It can measure between 0 – 9.999 % sulfur content in a sample. Consequently the equipment is applicable in more specific fields. Either a very small sulfur content as in diesel oil or a high sulfur content as in heavy oil can be determined.

■ Design

The SLFA-2800/2100 shows the new design philosophy of HORIBA. Also the function keys have been rearranged. Now the keys are arranged like cell phone keys which makes the operation and the data input even easier.



Easy to operate function keys.



Measurement results

Concentration of Sulfur in Samples (ppm)	Average Value(ppm)	Standard Deviation(σ) (ppm)
0	5.1	1.1
200	198.3	1.1
500	500.6	1.5

Sample: light oil Measurement conditions: measurement time 600 sec. Ten measurements



Features to ensure exact and reliable

■ Precision and repeatability

The new SLFA-2800/2100 investigates measurement results with an amazing precision. Due to the established miniature X-ray tube, which makes the purging of helium or another gas redundant extremely precise results can be achieved. In addition sudden fluctuations in temperature or detector time-drift are prevented through an automatic pulse-height compensation function. This system does also include a correction function for temperature and atmospheric pressure eliminates need for helium purging. As a result the measurement of low-concentration samples is no problem. The repeatability for the samples with a low sulfur content (0 %S) shows good performance with 1.6 ppm.

Also conce

■ A
Probl
are s
meas
of the

■ S
The
mode
cent
state

Lightweight and ultra-compact

Ultra-compact: 250(W) × 407(D) × 138(H)mm
Lightweight: 8kg

This small unit provides all of the necessary functions, even printing. The unit is easy to carry and ideally suited for on-site analysis.

High-precision measurements

Repeat accuracy: 1.5 ppm
Lower detection limit: 20 ppm

■ The SLFA series provides high precision and multiple functions

The SLFA Series is equipped with a multi-channel spectral analyzer for high-precision analysis and functions for multi-dimensional data processing. The SLFA Series can also perform calibrations for device drift caused by temperature fluctuations and atmospheric pressure fluctuations. Furthermore, there is no need to perform a helium purge.






■ Detector protection mechanism

The SLFA-2800/2100 includes a mylar film in order to protect the detector. Its design corresponds to the cell window's design and can be replaced easily. It is inserted between the detector and the sample cell. For that reason the detector is protected from sample drops in any case – even when sample cell breaks.

Cross section of sample cell and detector



Features for operator's convenience and safety

■ Easy-to-read large screen

Due to the large backlit LCD of the display the results of the measurements can be read well. Operations can be also easily monitored.

■ Coherent large printouts

The printouts are understandable and the size (80 mm wide) enables modest reading. All important information such as for example the sample ID, X-ray spectrum, date and time of the measurement, measuring time, number of accomplished measurements, average value, standard deviation and calibration curve graph are shown on the printouts.

```

**MEASURE
DATE: '01-07-30 14:27
CAL# : 1
MEAS TIME: 120 S
REPEAT : 3 TIMES
-----
SEQ# 1 POS#1
ID#HEAVY OIL1
1/ 3 10341.5 PPM
2/ 3 10338.0 PPM
3/ 3 10314.2 PPM
AVERAGE 10331.2 PPM
(K 0.51519)
STD DEV 14.3 PPM
(K 0.00053)
-----
SEQ# 2 POS#2
ID#LIGHT OIL1
1/ 3 98.0 PPM
2/ 3 101.1 PPM
3/ 3 103.1 PPM
AVERAGE 100.7 PPM
(K 0.10923)


```

Measurement result

```

**CALIBRATION-MAN2
DATE: '01-07-30 15:58
CAL# : 1
AUTO
QUADRATIC
% =A*K^2+B*K+C
B= 1.24783
C= -0.76994
D= 0.93330
RANGE: 1.06000 - 8.01990%
STD DEV 0.07522%
-----
STD# % K
01 1.06000 0.51449
(K 0.0003%)
02 4.26000 1.50892
(K -0.03614)
03 6.00220 1.83758
(K 0.11420)
04 8.01990 2.09990
(K -0.09181)

```



Auto-calibration

```

**SPECTRUM
DATE: '01-07-30 14:03
MEAS TIME: 10 S
-----
2000cps S
8
0 2 4 6 8 10
keV
S INT 29166:0 CFS
K(CORR) 0.51491
CORR COEF 1.00072
T 1603.5
S 13640.7 PPM

```

Spectrum

Features to ensure exact and reliable data

■ Precision and repeatability

The new SLFA-2800/2100 investigates measurement results with an amazing precision. Due to the established miniature X-ray tube, which makes the purging of helium or another gas redundant extremely precise results can be achieved. In addition sudden fluctuations in temperature or detector time-drift are prevented through an automatic pulse-height compensation function. This system does also include a correction function for temperature and atmospheric pressure eliminates need for helium purging. As a result the measurement of low-concentration samples is no problem. The repeatability for the samples with a low sulfur content (0 % S) shows good performance with 1.6 ppm.

Also measurements of samples with heavy sulfur concentration can be undertaken fast.

■ Automatic compensation of C/H ratio

Problems caused by measuring different oil types are solved by this function. Thus highly accurate measurement results can be achieved regardless of the oil type.

■ Spectrum measurement

The analyzer includes a spectrum measurement mode for assessments of the whole sample. In this mode you can print out a spectrum of the fluorescent X-rays, which is also useful for checking the state of the X-ray source or detector assembly.

■ Digital data output (RS-232C interface)

If storing of many results or further statistical processing of results is wished the measurement results can be put out to an external computer. Thus the analyzer can be used in laboratory or factory automation.


■ Safety mechanism

The operating and maintenance of our sulfur-in-oil analyzers is without any risks because no radioactive source is used. The exposure of X-rays is not possible due to a two-stage-fail-safe mechanism.

■ Turntable for sequential analysis and priority measurements (SLFA-2800)

The SLFA-2800 is equipped with a turntable, which makes sequential measuring of up to eight samples possible. With this feature the effort concerning each analysis decreases significantly. With the SLFA-2800 it is also possible to do priority measurements at any time. It does not matter if the SLFA-2800 is analyzing a single sample or a group of samples, it can always be paused to do measurements of higher priority.

■ System Diagram



A/DC : A/D converter
MCA : Multi-channel analyzer
AGC : Auto gain controller